

User Manual, English

EUROFLEX[®]

Sigma Work Chair



Eurovema

Contents	
INTRODUCTION	4
CONTACT DETAILS	4
USER INFORMATION	4
SERVICE AND WARRANTY.....	4
ASSEMBLY UPON DELIVERY	5
SAFETY RULES	5
SYMBOL KEY	6
WARNING	7
REPORTING OF ACCIDENTS AND INCIDENTS.....	8
CLEANING	8
STORAGE.....	8
MAIN LABEL.....	9
SECONDARY LABEL	9
CE MARKING	9
PRODUCT MARKING	10
CHAIR OVERVIEW	11
CHASSIS OVERVIEW	12
BACK HEIGHT	13
Side support, LS Back.....	13
BACK TILT, MANUAL BACK	14
BACK TILT - ELECTRONIC.....	15
BACK SUPPORT - SEAT DEPTH	16
SEAT ANGLE	17
SEAT TILT	18
ARMREST - SETTINGS	19
ARMREST - FIXED ANGLE.....	19
FOLDING ARMRESTS - SETTINGS	20
SEAT HEIGHT - MANUAL.....	21
Seat height in combination with seat rotation.....	21
SEAT HEIGHT - ELECTRIC.....	22
SEAT ROTATION.....	23
FOOTREST - HEIGHT.....	23
LEG SUPPORT - FULL FOOTPLATE TILT	25
FOOTPLATE - TILT	26
BIFURCATED FOOTREST - HEIGHT ADJUSTMENT	27
BIFURCATED LEG SUPPORT - TILT.....	28

BIFURCATED FOOTPLATE - FOOT TILT	29
FOOT RING.....	30
BRAKE - CHASSIS-MOUNTED	31
CHARGING - GENERAL	32
CHARGER	Fel! Bokmärket är inte definierat.
BATTERY LOCATION.....	33
BATTERY FUNCTION	34
REPLACING THE BATTERY	34
TECHNICAL DATA - SIGMA.....	35
ENVIRONMENT	36
TECHNICAL DATA - CHASSIS	36
TECHNICAL DATA - ELECTRONICS.....	36
TECHNICAL DATA - SEAT, BACK SUPPORT & ARMRESTS.....	37
ACCESSORIES LIST	38
DESTRUCTION.....	39

INTRODUCTION

We congratulate you on your choice of work chair! We hope you will be pleased with this Euroflex product from Eurovema, designed and built in Gislaved, Sweden. The work chair is designed to satisfy very high demands for ergonomics, sitting comfort and function. The chair has many options for adjustment enabling it to be adapted to every need. Read the user manual carefully so you can make use of all the options provided by your chair. Eurovema reserves the right to make changes to this manual and its contents.

CONTACT DETAILS

Manufacturer



Eurovema Mobility AB
Baldersvägen 38
SE-332 35 Gislaved, Sweden
Telephone +46 371 390 100
<http://www.eurovema.se>
info@eurovema.se

USER INFORMATION

- This work chair is classified as a class 1 medical device.
- The product is intended to be used by people who require assistance to move or carry out activities whilst seated and need support to get out of a chair.
- The product is only intended for indoor use.
- There are no known contraindications.
- Maximum user weight is 80 - 160 kg, depending on the seat system used.
- To adapt the chair to heavier users, please contact Eurovema.
- When used in accordance with this User Manual, the expected service life of the product is ten years.
- The product has been developed and tested by Eurovema, and subjected to third-party testing by an independent testing body. Eurovema is responsible for ensuring that the product is safe to use.



This symbol warns of situations that may pose risks to the user.



Read the User Manual carefully.


SERVICE AND WARRANTY

If the work chair is used on a daily basis, it should be sent to for inspection once every other year. This is to check that the function and safety of the chair are maintained during its entire service life. Expected service life is 10 years if used in accordance with our instructions, stipulated maintenance intervals, and intended use. If your chair requires servicing, please contact your local Technical Aid Centre.

We give a 4-year warranty on the metal structures in our products, a 2-year warranty on other parts except upholstery, wheels, and batteries, which come with a 1-year warranty, unless otherwise agreed. Normal wear and tear is not covered by the warranty. We recommend that our customers to use the product in accordance with the user manual. For instructions in respect of servicing and reconditioning, **see Service Manual.**

ASSEMBLY UPON DELIVERY

- Open the packaging and check that no damage has occurred during transit. The work chair is supplied with the back support and folding armrests not fitted.
- The armrests are equipped with stop screws, and the back support features a stop ball that must be pressed in during assembly and disassembly to prevent them from releasing when adjusting the back rest to its end position.
- Fit the back support and armrests as shown in the figures and tighten the lock knob (1).
- Fit stop screws (2) using a Phillips screwdriver.





































SAFETY RULES

- Read the User Manual carefully before using the chair.
- Perform a visual inspection of the chair every day before use.
- The work chair is intended to be used on a level surface in a normal indoor climate.
- Take care when adjusting the manual seat tilt while sitting in the chair as there is a risk you could fall out of the chair.
- If you find that the chair has suffered damage or identify any changes in its function, contact the service organisation (Technical Aid Centre) immediately.
- If the chair is to be used on a regular basis, it should be sent routinely to the service organisation for inspection.
- Make sure you tighten the lock knob and screws properly after making adjustments.
- Pay attention to screws or parts that become loose or fit loosely as they may impact safety. Contact the service organisation immediately.
- Service and maintenance must only be carried out by a Eurovema authorised technician.
- Only original parts from Eurovema may be used.
- The chair must not be equipped with any accessories or components other than those approved by Eurovema.
- Do not exceed the maximum stated user weight.
- Check that stop screws are fitted to the armrests and a stop bearing in the back support attachment, this prevents them from coming loose when adjusting the back support to its end position.
- Users may use the functions mentioned below in this manual. Users may not use functions not mentioned in this manual.
- Under no circumstances must the user sit in the chair whilst it is being serviced, except when armrests and leg supports are being replaced and steps can be taken to ensure that the user does not fall out of the chair.



Chairs with adjustable seat tilt must be in the fixed position and the brake applied when the user gets in or out of the chair.

SYMBOL KEY

	Warning		Read User Manual
	Warning, risk of crushing		Seat tilt
	Do not iron		Adjustable back tilt
	Do not use dry cleaning solvents stronger than tetrachloroethylene		Back height
	Do not use bleach		Seat depth
	Do not tumble dry		Armrest height
	Machine wash 60		Operation button, electric function
	Protect against liquid		Charging port
	Handle with care		Input control
	This way up		Safe Working Load
	Permitted temperature		Applied Part Type B
	Permitted relative humidity		
	Permitted height above sea level		
	Manufacturing date		
	Manufacturer		
	For indoor use only		
	Maximum user weight		
	Serial number		
	CE mark shows Compliance with European regulations		
	Electrical components are to be sent for dedicated recycling		
	The product is part of a recycling system		
	Component weight exceeds 10kg		
	Operation cycle		

WARNING



A product that exhibits impaired or changed performance must be taken out of use immediately. Contact the service organisation (Technical Aid Centre) or Eurovema AB immediately. The product may not be used again until it has been inspected by an authorised technician.



The product may not be subjected to a load in excess of the stated maximum user weight as this may result in components being damaged, which may lead to accidents.



Check that the parking brake is activated when getting in or out of the chair. If it is not, there is a risk that the product will roll unintentionally, which may result in you falling.



Modifications to the product may only be made by an authorised technician in consultation with the prescriber. Should a custom adaptation be required, contact a Technical Aid Centre or Eurovema Mobility AB for advice. Modifications to the product will invalidate the CE mark and Eurovema will no longer accept any liability. Modifications carried out by an authorised technician that have been agreed with Eurovema may allow the CE mark to remain valid.



Repairs and technical service may only be carried out by personnel authorised by Eurovema Mobility AB in order for the CE mark to remain valid. Repairs and technical maintenance may not be carried out with the user sitting in the chair. For more detailed servicing and maintenance instructions, see the Service Manual.



The product is intended to be used by the person and for the purpose for which it has been prescribed. The product has been configured for this person and its settings may not be changed by anyone other than the person who prescribed the product. If the product is to be transferred to another user, the chair must undergo a new test.



In certain circumstances, a product with an electric seat lift may elevate the seat to a height that will result in the user being unable to reach the floor with their feet, thereby creating the risk of a fall should the user attempt to get out of the chair from this position. Make sure that personnel or another person is present if the highest seat elevation setting is used.



Components and parts must be fitted in accordance with the assembly instructions in order to ensure that they are secure and do not come loose. All cables must be secured using cable ties.



In order to maintain the validity of the CE mark, the product may not be equipped with any components or accessories other than those approved by Eurovema Mobility AB.



Only electrical components specifically intended for use with the product may be connected to the product's control system. Components must be connected to the control system by an authorised technician.



In order to maintain the validity of the CE mark, and to ensure compliance with the

terms of the warranty, only chargers and batteries specifically designated for the product and supplied by Eurovema Mobility AB may be used.



A product with electrical functions may experience interference from mobile telephones or other RF communication equipment.



The product is features many small parts and screws that, if they come loose, may pose a choking hazard for children and pets.



In certain configurations, the product is equipped with cables and wires that functional reasons, cannot be completely secured and, as a result, may constitute a risk of strangulation.



On products with an electric brake, it is important the the battery is not completely discharged as this will result in the brake ceasing to work. In such situations, and if the brake is not already applied, the chair may roll when you sit down or get up, with a risk of falling as a result.

REPORTING OF ACCIDENTS AND INCIDENTS

EU-based manufacturers of medical devices are obliged to systematically track how their products work in practical use.

The awarding of a CE mark shows that the product has undergone a series of risk analyses and tests, with steps taken to minimise any risks identified as far as possible. If, despite everything, accidents or incidents occur, such events must be reported to Eurovema Mobility AB and the relevant national authority.

CLEANING

Work chair surfaces are to be cleaned using a damp cloth used in conjunction with a mild detergent, e.g. washing-up liquid. Place the chair on a stable surface and apply the brake. Seat upholstery can be cleaned using commercially available upholstery cleaners. Loose covers can be machine washed at 60°.



NB Cleaning must be done in a well ventilated are and not in the presence of any naked flames.


STORAGE

The work chair must be stored in a dry environment at room temperature. Keep the chair away from high temperatures, intense cold, and strong sunlight. Metal surfaces can become very hot if exposed to sunlight. In instances of intense cold, allow the chair to come to room temperature before use. In addition, do not expose the chair to water, other liquids, or chemicals.

MAIN LABEL


The label on the product contains important information.

- 1) Product name
- 2) Part number
- 3) Serial number UDI-PI
- 4) Barcode
- 5) Max. user weight
- 6) Medical device
- 7) CE marking
- 8) Read the manual before use
- 9) Manufacturing date
- 10) Manufacturer's name
- 11) Applied Part Type B



SECONDARY LABEL

- 1) Applied part type B
- 2) Intermittent function 2 / 18 min
- 3) Permitted room temperature
- 4) Permitted relative humidity
- 5) For indoor use only
- 6) Permitted height above sea level
- 7) Safe Working Load



CE MARKING

The product is CE-marked in accordance with the Medical Devices Regulation (MDR 2017/745) of the European Parliament. The CE mark can be found on the label. The product is compliant with the requirements of the standard.

EN1335-2:2018 Office furniture. Office work chair-Safety requirements

EN 60601-1:2006 Medical electrical equipment - part 1: General requirements for basic safety and essential performance.

EN 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances.

SS-EN 1021-2:2014 Furniture - Assessment of the ignitability of upholstered furniture - Part 2: Ignition source: Gas flame equivalent to a lit match.

SS-EN 1021-1:2014 Furniture - Assessment of the ignitability of upholstered furniture - Part 1: Ignition source: Burning cigarette.

SS-EN ISO 21856:2022 Assistive products – General requirements and test methods (ISO21856-2022)

SS-EN ISO 14971:2020 Application of risk management for medical devices.

PRODUCT MARKING

Loose covers can be machine washed.

INTERMITTENS
2 / 18 MIN
DUT CYCLE





Crush risks are reported in connection with instructions for each function that may constitute a risk of crushing. Marking appears on the chair where there is a risk of crushing.

CHAIR OVERVIEW

- 1) Seat
 - a) SitRite
 - b) Forma
- 2) Back support
 - a) SitRite
 - b) SitRite LS
 - c) Forma
- 3) Armrests
- 4) Chassis
- 5) Brake lever
- 6) Gas Seat Lift
- 7) Electric Seat Lift
- 8) Foot ring




|



CHASSIS OVERVIEW

- 1) Front wheel
- 2) Foot ring
- 3) Parkring brake lever
- 4) Chassis
- 5) Brake rod
- 6) Rear wheel
- 7) Lifting pillar, Electric
- 8) Seat lift attachment, Gas




BACK HEIGHT

The back support is height adjustable using a screw and a lever.

- 1) Loosen the lever a ½ turn by turning it anticlockwise.
- 2) Set desired height.
- 3) Lock the back by turning it clockwise a ½ turn.
- 4) To adjust the height further, loosen the adjustment screw a few turns.
- 5) Set desired height.
- 6) Lock adjustment screw.


The screw can be replaced with a locking wheel.



Side support, LS Back

The side supports are adjustable horizontally.

- 1) Loosen the screws (1) a few turns using a 4mm Allen key.
- 2) Adjust to the desired width.
- 3) Tighten the screws.



BACK TILT, MANUAL BACK

Back tilt

The back tilt and back cushion are individually adjustable.


- 1) Undo the safety screw and lever to tilt the entire back mechanism backwards or forwards.
- 2) Set desired back angle.
- 3) Lock safety screw and lever.

Back cushion

- 4) Loosen the lever by turning it anticlockwise a ½ turn.
- 5) Set the desired angle of the back cushion.
- 6) Lock the lever by turning it clockwise a ½ turn.



Adjustment may pose a risk of crushing



BACK TILT - ELECTRONIC




The back mechanism tilt is electronically adjustable.

- 1) Tilt back by moving the rocker switch upwards.
- 2) Tilt forwards by moving the rocker switch downwards.

An optional hand control is also available, functions as shown in the figure below.
The back tilt can be limited electronically.



Do not adjust the back tilt in the immediate vicinity of another person.



BACK SUPPORT - SEAT DEPTH

Seat depth can be adjusted using the back support. Set the desired seat depth by moving the back support backwards or forwards. Fine tuning of the seat depth can also be done using the back tilt.

- 1) Loosen the locking screw/knob under the seat on the rear edge.
- 2) Push in or pull out the back mechanism to the desired seat depth.
- 3) Tighten the locking lever/screw.
- 4) Fine tune using the back tilt, see previous section *Back Tilt*



SEAT ANGLE

The seat angle is adjustable, and several different tilt mechanisms can be used:


1. Turnbuckle
2. Crank
3. Gas piston
4. Electric tilt.

1) The **turnbuckle** is adjusted using a 17mm box wrench.

2) **Crank tilt** is adjusted by turning the wheel.



Adjusting the seat tilt may pose a risk of crushing, so do not carry out this procedure with children or pets in the immediate vicinity.



SEAT TILT

3) Gas-controlled seat tilt is regulated using the control located under the armrest cushion. The control can be moved from right to left side. To adjust the seat angle, press the control upwards towards the armrest cushion and release once the desired angle has been reached.

4) Electric seat tilt is controlled using either the rocker switch located under the armrest or the hand control (see sections *Back Tilt - Electric*, *Seat Height - Electric* for more details). Electric seat tilt can be limited electronically.



Adjusting the seat tilt may pose a risk of crushing, so do not carry out this procedure with children or pets in the immediate vicinity.




ARMREST - SETTINGS

- 1) To adjust the height of the armrests, loosen the wheel (1) a few turns clockwise and set to the desired height.
- 2) Tighten the wheel.
- 3) To adjust armrest width, loosen the lever (2) a ½ turn clockwise and adjust the the desired width.
- 4) Tighten the lever.



ARMREST - FIXED ANGLE

The armrests can be fitted in a position that does not change when the seat is tilted if the seat is supplemented with attachments for fixed angle armrests. Set width and height as shown in the figure above.




FOLDING ARMRESTS - SETTINGS

The armrest can be folded backwards by pulling the spring-loaded control (1) upwards and pushing the armrest backwards at the same time.

- 1) The height of the armrest is adjusted by the lever or screw (2) being loosened a $\frac{1}{2}$ turn. If the armrest is equipped with a screw, use Allen key no. 5.
- 2) Tighten the lever or screw once the desired height has been reached.
- 3) For width adjustment, see the previous section.



Make sure that no-one can be crushed when the armrest is lowered.



SEAT HEIGHT - MANUAL

- 1) The seat is raised or lowered by moving the control under the armrest (2) upwards. To lower the seat, sit in the chair and keep the control pulled up. The seat will drop using your body weight.
- 2) Release the control once the desired height has been reached.




When raising the seat with the chair under a table/desk there is a risk that your hand or legs may be crushed.



Seat height in combination with seat rotation

On chairs equipped with seat rotation, the seat height control is located under the seat, see (1)

- 1) Lift the lever to activate the seat lift.
- 2) Release the the lever once the seat has reached the desired height.




SEAT HEIGHT - ELECTRIC

Seat height is controlled electronically using the rocker switch located under the armrest. Pull the rocker upwards to raise the seat and press it down to lower the seat.

A hand control is available as an optional extra with functions as shown in the figure below.




Take care when lowering the seat to ensure that there is no-one gets crushed. When raising the seat with the chair under a table/desk, there is a risk that your hands or legs may be crushed.



SEAT ROTATION

Seat rotation is available as an optional extra and allows the seat to rotate 90° in both directions. To rotate the chair, apply the brake and move the control upwards at the same time you turn the seat. The control can be released once you have started to turn and will fall back to the starting position once the seat has turned 90°.




FOOTREST - HEIGHT

- 1) Loosen the safety screw (1) using Allen key no. 5
- 2) Loosen the screw (2) a few turns and adjust to the desired height.
- 3) Tighten the screw (2) and the safety screw (1).



Do not stand on the footplate, the chair will tip forwards




LEG SUPPORT - FULL FOOTPLATE TILT

- 1) Loosen the screw (1) a few turns using Allen key no. 5.
- 2) Set desired angle.
- 3) Tighten the screw.



Do not stand on the footplate, the chair will tip forwards




FOOTPLATE - TILT

The footplate can be adjusted horizontally. To tilt it downwards, screw the screws anticlockwise using Allen key no. 5 until the desired angle is reached. Lock the screws using nuts. To tilt upwards, undo the nuts and unscrew anticlockwise.



Do not stand on the footplate, the chair will tip forwards



BIFURCATED FOOTREST - HEIGHT ADJUSTMENT

- 1) Loosen safety screw (1) and screw (2) a few turns so that the footplate can be raised or lowered.
- 2) Once the desired height has been reached, tighten screws (1) and (2).



Do not stand on the footplate, the chair will tip forwards




BIFURCATED LEG SUPPORT - TILT

The leg supports can be angled. Loosen the screw a few turns using Allen key no. 5, adjust to the desired angle, and tighten the screw.



Do not stand on the footplate, the chair will tip forwards




BIFURCATED FOOTPLATE - FOOT TILT

The angle of the footplate can be adjusted to enable the user to find a comfortable position. To do this, loosen the screw a few turns using Allen key no. 5 and adjust the angle of the footplate (1). Tighten the screw. The footplate can be lifted without tools (2) to facilitate getting in and out of the chair.



Do not stand on the footplate, the chair will tip forwards



FOOT RING

The foot ring can be folded backwards, resting against the chassis in the fully folded position.



Fold the foot ring away before getting out of the chair in order to reduce the likelihood of a fall.



BRAKE - CHASSIS-MOUNTED

Apply the brake by moving the brake lever forwards, and release the brake by moving the lever back again.



Always apply the brake before getting in or out of the chair.



CHARGING - GENERAL

To ensure full battery performance for as long as possible, it is important that the battery is charged regularly. Usually, it is sufficient to charge the battery 1–2 times a week, preferably during the daytime when the chair is not in use. Do not interrupt the charging process until it has completed. Charging is complete once the LED on the rocker switch box shows a constant orange. When the chair is in long term storage, the battery should be charged once every four weeks. During charging, place the chair so that the charging cable can be easily connected and disconnected. Whilst the chair is charging, only the electric brake will be operational (engage/disengage) in order to allow the chair to be moved without disconnecting the charger, no other function is usable.



Do not leave the chair charging for long periods of time. Do not leave the charger connected to the wall outlet once charging is complete.




On chairs with electric functions, it is particularly important to charge the battery on a regular basis in order to ensure that the electric brake can be applied and to control seat position so that the user can get in and out of the chair safely.

CHARGER

TP7M

The battery must be charged when the battery issues an audible signal and the LED (4) on the rocker switch box flashes when activating any of the chair's electronic functions. In order to protect the battery, charging must be done as soon as this audible signal is heard and the LED starts to flash.

- 1) Connect the charging plug (2) to the socket (1) in the rocker switch box.
- 2) Connect the charger wall plug to a 220V wall outlet. The charger LED (3) will be green. The LED (4) will flash for a few seconds and then go out.
- 3) Charging is complete once the orange LED (4) on the rocker switch box is a constant orange.
- 4) Disconnect the charger from the wall outlet, and then disconnect the charging plug from the chair.



BATTERY LOCATION

The battery **on/off** buttons and **LEDs** can be found on the top side of the battery box, see figure.




BATTERY FUNCTION

When the charger is not connected

Button	Function	Activity
On	Press and hold for more than 1 second	Battery activated
Off	Press and hold for more than 3 seconds	Battery deactivated

Status	Capacity	LED	
Charged	>80%	Green	Constant light
Loading	<80%	Green	Rapid flash, 1/1 second
Charging required	>25%	Green	Long flash, 0.5/4 seconds
Battery low	<25%	Orange	Long flash, 0.5/4 seconds
Protection mode	<20%	-	Off
Fault warning	-	Orange	Constant light

To disengage electrical functions, press the OFF button located on the battery. To reinitialise functions, connect the cable to the charging port in the seat until you hear an audible signal. For more detailed instructions, see page 32



See technical data

REPLACING THE BATTERY

To replace the battery, the entire battery box must be replaced. For instructions on how to do this, see the installation procedure outlined in Replacing the Battery.



Repairs and other procedures on the battery may only be carried out with the approval of Eurovema. The battery box contains no user serviceable parts and must not be opened. Only Eurovema-approved batteries may be used in the product.

TECHNICAL DATA - SIGMA

Class	Work chair ISO 180903
Weight, electric	42kg
Weight, manual	34kg
Weight SWL, total weight max user weight + chair	175kg
Maximum user weight	160kg
Battery, electric version	TBB4 24V 2Ah integrated monitoring for more info. See Data TBB4

Seat height EL with single wheel	Without rotation	With rotation
Seat height EL 100mm, clearance 50mm	383-483mm	393-493mm
Seat height EL 100mm, clearance 70mm	403-503mm	413-513mm
Seat height EL 180mm, clearance 50mm	463-643mm	473-643mm
Seat height EL 180mm, clearance 70mm	483-663mm	493-673mm
Seat height EL 250mm, clearance 50mm	533-783mm	543-783mm
Seat height EL 250mm, clearance 70mm	553-803mm	563-813mm

Seat height EL with double wheel	Without rotation	With rotation
Seat height EL 100mm, clearance 50mm	363-463mm	373-473mm
Seat height EL 100mm, clearance 70mm	383-483mm	384-493mm
Seat height EL 180mm, clearance 50mm	443-623mm	453-633mm
Seat height EL 180mm, clearance 70mm	463-643mm	473-653mm
Seat height EL 250mm, clearance 50mm	513-763mm	523-773mm
Seat height EL 250mm, clearance 70mm	533-783mm	543-793mm

Seat height, Gas, regardless of wheel	Without rotation	With rotation
Seat height, Gas 95mm, clearance 30mm	311.5-406.5mm	321.5-407.5mm
Seat height, Gas 155mm, clearance 30mm	408.5-563.5mm	418.5-573.5mm
Seat height, Gas 200mm, clearance 30mm	468.5-668.5mm	478.5-678.5mm

Seat tilt	Front	Rear
Seat tilt electric, crank, turnbuckle, and gas	-14°	+16°

Back angle	Front	Rear
Seat tilt, electric, manual, and gas	-23°	+7°

ENVIRONMENT

Condition	Temperature	Humidity	Atmospheric pressure
Normal use	5°C - 40°C	10% - 93% RH	50kPa – 106kpa
Storage and transport	-5°C - 40°C	10% - 93% RH	50kPa – 106kpa

TECHNICAL DATA - CHASSIS

Data	Facts
Width, chassis with single wheel	582mm
Width, chassis with double wheel	603mm
Length, chassis with single wheel	592mm
Width, chassis with double wheel	603mm
Chassis weight	9.9kg
Manual brake	Chassis-mounted lever brake
Electric brake	Timotion TL16 500N stroke 50mm
Wheel diameter	100mm

TECHNICAL DATA - ELECTRONICS

Data	Facts
Input voltage	30V DC Max 6A
Output voltage	30V DC 6A
Work cycle	10% max 2min on / 18min off
Voltage	24V DC
EMC	Class 2
Protective enclosure	IP66
Electric protection class	Class 2
Battery TBB4	24V 2Ah AGM



Chairs that are equipped with electronic components may be disrupted by radiowaves, e.g. mobile telephones.



Chairs equipped with electronic components are classed as Applied Part Type B, protection class II with internal power source.

TECHNICAL DATA - SEAT, BACK SUPPORT & ARMRESTS

Data	Facts
Seat (width x depth)	40x40, 40x46, 45x46, 45x50, 50x46, 50x50cm
Seat width (between armrests)	38-56cm
Seat depth	40-59cm
Back support, SitRite (height x width)	29x32, 36x40, 42x45, 47x47cm
Back support, SitRite LS	17x37, 24x43, 28x45, 33x47cm
Back support, Forma	37x23, 37x46, 37x48cm
Back support, Classic	35x30, 35x44, 40x30, 40x44cm
Adjustable back tilt	30°
Armrest, height adjustable (follows tilt)	0-26cm
Armrest, height adjustable (does not follow tilt)	^a 0-26cm (Theoretical -15.5-26cm)
Armrest plates	25x8x3, 30x8x3, 40x8x3, 25x8x5, 30x8x5, 40x8x5cm
Seat, Coxit (width x depth)	45x45, 50x50cm
Seat rotation	180° lockable in five positions
Upholstery	Black Atlantic

a) Requires long armrest post.

ACCESSORIES LIST

The following accessories can be replaced in order to change the functionality of the product:

1. Control rocker switch, available in hand control format, wired and wireless.
2. Wheels, available in single and double wheel variants.
3. Brake, available as manual brake and electronic brake.
4. Lifting pillar, available in 3 different heights.

Contact Eurovema for a complete list of parts that can be repaired/replaced - authorised service personnel only.

DESTRUCTION

Products that are to be permanently put out of use must be disassembled and sorted at source. Assembly and disassembly instructions can be found on the website in the *Assembly Instructions* section. Detailed diagrams can also be found on the website in the *Exploded Diagrams* section. Point your browser to www.eurovema.se

Private individuals should contact Eurovema for more information. Eurovema will assume responsibility for destruction.

Disassembled components are to be sorted in accordance with the **table below**.



The disassembly process is not risk-free. There is a chance you may suffer crush injury or burn injury if you do not follow the procedure correctly.. Disassembly must be carried out by an authorised technician. P



Metal	Electronics	Combustible
Chassis	Control electronics	Seats
Seat cross	Connection boxes	Back support
Back mechanism	Rocker control	Armrest cushions
Armrest posts	Battery box (batteries must be isolated)	Neck support cushions
Leg support	Electric actuator	Detachable upholstery covers
Accessory attachments	Lifting pillar, Electric	Trunk support cushions
Seat lift, Gas	Cables	Other filling materials
Foot ring		

Batteries

Discarded lead batteries must be sent to a recycling centre.

