

i4 SitRite

Euroflex i4 is the smooth, 4-wheel electric indoor wheelchair that is the much-awaited successor to the popular Miniflex and is compatible with the SitRite seat system. With a lifting height of 30 cm, which affords a maximum seat height of 73 cm, the i4 SitRite allows users to get in close to tables and work surfaces. The forward tilt makes it a chair that is easy to get in to and out of, whilst also facilitating activities near a work surface. i4 is a front wheel drive chair, and the slightly larger drive wheels makes crossing thresholds a breeze. The articulated rear wheels allow the frame to flex, keeping the chair in contact with the floor.

Type of impairment

- minor
 moderate
 extensive
 indoor
 outdoor
 child
 adult



SEAT SYSTEM

SitRite / SitRite LS



The inherent benefits of SitRite are concealed within the seat, taking the form of two foam materials of different densities. The lower, harder foam section helps to hold the pelvis in a good position, while the upper, softer foam section distributes pressure by being thicker under the ischium than it is under the thighs.

SitRite's back is ergonomically designed, ensuring that a large part of the trunk is supported without impairing arm freedom of movement. The SitRite LS (Lateral Support) variant features a more enclosed back. SitRite LS has been produced specifically for people with trunk muscle weakness and correctable scoliosis. The enclosing sides are structured around polyester covered plates that can be adjusted whilst the user is in the sitting position.

Both SitRite and SitRite LS can be fitted with a lumbar pillow to provide extra support.

UPHOLSTERY



Black fabric

Seat (width x depth)	36x40, 36x45, 40x40, 40x45, 45x45, 45x50, 50x50 cm
Seat width (between armrests)	36-58 cm
Seat depth	38-56 cm
Seat height (height to underneath the seat)	43-73 cm
Seat tilt, electrical, crank & turnbuckle	-15 to 24
Back support, height, adjustable	37-62 cm
Back (width x height) - SitRite - SitRite LS *	35x43, 40x45, 45x47 cm 24x43, 28x45, 33x47 cm
Adjustable back angle	30°
Armrest, height, adjustable	0-26 cm
Armrest platform (length x width x thickness)	30x8x3, 30x8x5, 40x8x3, 40x8x5 cm

*) The width of SitRite LS is measured inside the narrowest part between the back rest wings. The side plates on the SitRite LS can be adjusted in width so that the specified width can be 3-5 cm wider depending on the size of the back rest.

CHASSIS

i4



The market's most compact four-wheel chair chassis for indoor use is no larger than a traditional work chair. The easily manoeuvred and compact design makes it easy to get close to tables or work surfaces. The chassis has front-wheel drive with a jointed rear axle, which ensures that the rear wheel pair is always in contact with the floor.

Chassis, width	56 cm
Chassis, length	73 cm
Total weight	75 kg
User weight, max.	130 kg
Turning circle	55 cm
Max speed	4 km/h
Max distance	approx. 15 km
Electronics	PG
Obstruction clearance	5 cm
Batteries	2x12 V/33 Ah