


Rollators made in Sweden



VOLARIS

Stable and easy to operate



Our rollators are not only developed in Sweden, at Eurovema Mobility we also make and assemble them. You can rest assured in the knowledge that each and every Volaris rollator is based on tried and tested technology that has a reputation for high quality, functionality, and user-friendliness.

Volaris rollators are easy to operate yet stable. The cross-shaped flexing chassis enables the rollator to move smoothly and stably, making it easy for the user to walk over uneven surfaces. All models can be easily collapsed and stowed vertically. Thanks to its low weight and collapsible design, Volaris is easy to transport anywhere by car.

All available adjustments, including handle height, are clearly marked with orange coloured controls, a feature that is highly appreciated by many users. Brake cables are integrated into the frame, eliminating the risk of them catching on anything. All Volaris rollators are equipped with a stable seat that is also height-adjustable. All models are also available in extra wide versions that accommodate users weighing up to 160 kg and feature wider seats.

» **VERY SAFE AND STABLE, BOTH WHEN USED TO WALK WITH, AS WELL AS WHEN SITTING DOWN AND GETTING UP.**



COMPACT

A uniquely smooth indoor rollator that still offers reliable stability. With a total width of only 53 cm, users can easily get around their home.

» SMALL AND SMOOTH INDOORS



SMART

A popular all-rounder that features medium-large wheels that allow users to go anywhere, both indoors and outdoors.

» ONE OF THE LIGHTEST ON THE MARKET



PATROL

A much appreciated and stable lightweight rollator with four large wheels that can cross even the most uneven outdoor surfaces with ease. The most adjustable handles on the market - up to a height of 107 cm. All Patrol models can handle user weights of up to 160 kg.

» THE PERFECT CHOICE FOR TALL PEOPLE



DISCOVERY

Equipped with softer and wider wheels that minimise handle vibration. Specifically developed for people who suffer from rheumatism or experience pain in the hands.

» ADAPTED TO MEET THE NEEDS OF RHEUMATISM SUFFERERS



VOLARIS RA

RA-support with underarm supports that are adjustable in all directions, providing safe and stable use. Also available with extra large wheels for even more stable use outdoors.



VOLARIS HP

Designed to act as a walking aid for people with weakness on one side of the body (Hemiplegia). The HP support can also be replaced with an underarm support similar to RA. For user weights up to 160 kg.



VOLARIS MINI

A small, 53 cm wide rollator for both indoor and outdoor use that is designed for children and shorter adults.



VOLARIS KID

Designed for children. Smaller wheels and a lower seat and handle height than the Mini.

Eurovema

Eurovema Mobility AB
Baldersvägen 38 | SE-332 35 Gislaved | Sweden
info@eurovema.se | eurovema.com
Phone: +46 371 390 100