

SitRite Seat System

The SitRite seat system is designed to provide users with additional support and comfort. It has been produced in collaboration with seat ergonomists and has proved to be very popular for many years. The range now also includes SitRite LS (Lateral Support), which provides a more enclosing back piece. SitRite LS has been produced specifically for people with trunk muscle weakness and correctable scoliosis. The enclosing sides of SitRite LS are built around polyester coated plates, which can be adjusted whilst the user is in sitting position.

Both SitRite and SitRite LS can be fitted with a lumbar pillow to provide extra support. The lumbar support is fitted between the cushion and back plate. The support is inflated using a hand pump.

The SitRite seat consists of two layers of foam material, each of which has a different density and performs a different function. The lower, harder layer helps maintain the pelvis in a good position. The upper, softer foam part helps distribute pressure by being thicker under the ischium and thinner under the thighs.

SitRite and SitRite LS can be used across the entire range of work and electric wheelchairs, i.e. Junior, Sigma, Flexbase, Hi/low, Miniflex, and Flexmobil. The SitRite LS back can also replace an existing SitRite back.

The range includes seat and back covers for all models and in all sizes.



SITRITE SEATS



Seat model	Part no.	Inner dimensions (W x D cm)	Intended use
SitRite 26x31	624570-2631	26x31	Junior
SitRite 29x32	624570-2732	27x33	Junior
SitRite 32x36	624570-3236	31x37	Junior
SitRite 36x40 Junior	624570-3640	34x41	Junior
SitRite 40x40 Junior	624570-4040	38x41	Junior/adult
SitRite 36x45	624570-3645	34x46	Junior/adult
SitRite 36x40	624520-3640	34x41	Adult
SitRite 40x40	624520-4040	38x41	Adult/Junior
SitRite 40x45	624520-4045	38x46	Adult/Junior
SitRite 45x45	624520-4545	43x46	Adult
SitRite 45x50	624520-4550	43x50	Adult
SitRite 50x50	624520-5050	48x50	Adult

The inner dimensions of the seat are measure in relation to where the seat is shaped, i.e. the seat indent. See marking on the seat for "inner dimensions".

SITRITE BACKS



Back model*	Part no.	Inner dimensions (W x D cm)	Intended use
SitRite 25x32	224576-2532	25x32	Junior
SitRite 27x37	224576-2732	27x37	Junior
SitRite 35x43	224526-3640	35x43	Junior/Adult
SitRite 40x45	224526-4245	40x45	Adult
SitRite 45x47	224526-4550	45x47	Adult
SitRite LS 14x32	224536-2532	14x32	Junior
SitRite LS 17x37	224536-2732	17x37	Junior
SitRite LS 24x43	224536-3640	24x43	Junior
SitRite LS 28x45	224536-4245	28x45	Junior/Adult
SitRite LS 33x47	224536-4550	33x47	Adult

* All models include a back profile.

The inner dimensions of the back on SitRite LS are measured at the narrowest part between the wings. See marking in the illustration.

The width of the side plates on SitRite LS can be adjusted, meaning that the specified width may be a couple of centimetres narrower or a couple of centimetres deeper.