



Important information for users of Flexbase SitRite / Flexbase Forma / Clinic work chair with chassis as shown in figure above

It has been discovered that chairs sold during the period 05.08.2019 to 31.12.2021 were fitted with leg halves that contained a manufacturing defect that may cause the screw holding the leg halves together to fail. In the event of this screw failing, there is significant risk of the user falling out of the work chair, thereby exposing themselves to physical danger.

Please notify all affected individuals in your organization.

We will be contacting you during November to discuss what will happen next.

The plan is to carry out works on affected chairs during 2025.

Important Safety Notice to the Market

Safety Notice FSN 104121
2024-11-18

If you have any questions, please contact Eurovema Customer Services, either by e-mail at info@eurovema.se or telephone +46 371 39 01 00

Quality Manager: peter.kjell@eurovema.se

Affected serial numbers:

1040251938	1040253998	1040256700	1040260027
1040251944	1040254011	1040256701	1040260196
1040252007	1040254031	1040256703	1040260200
1040252481	1040254062	1040256705	1040260391
1040252544	1040254239	1040256707	1040261015
1040252548	1040254278	1040256709	1040261141
1040252760	1040254281	1040256766	1040261275
1040252783	1040254400	1040256769	1040261287
1040252955	1040254401	1040256772	1040261583
1040252956	1040254402	1040256884	1040261588
1040252958	1040254605	1040256888	1040261592
1040253155	1040254668	1040257148	1040261619
1040253156	1040254669	1040257217	1040261951
1040253160	1040254728	1040257254	1040261957
1040253203	1040254729	1040257255	1040261965
1040253208	1040254730	1040257353	1040261967
1040253411	1040255031	1040257420	1040262152
1040253412	1040255076	1040257421	1040262158
1040253413	1040255245	1040257718	1040262237
1040253414	1040255254	1040257719	1040262289
1040253421	1040255287	1040257782	1040262521
1040253424	1040255339	1040257783	1040262526
1040253425	1040255374	1040257841	1040262640
1040253436	1040255377	1040257842	1040262830
1040253440	1040255596	1040257853	1040262959
1040253441	1040255703	1040257876	1040263167
1040253548	1040255728	1040257877	1040263168
1040253638	1040255779	1040258080	1040263169
1040253639	1040255781	1040258133	1040263170
1040253641	1040255818	1040258137	1040263375
1040253642	1040256015	1040259525	1040265506
1040253643	1040256019	1040259546	1040265772
1040253677	1040256020	1040259565	1040265826
1040253730	1040256046	1040259694	1040265830
1040253804	1040256047	1040259749	1040265889
1040253805	1040256367	1040259755	1040265890
1040253806	1040256394	1040259760	1040266151
1040253986	1040256439	1040259852	1040266203
1040253997	1040256440	1040260019	1040266222

Important Safety Notice to the Market

Safety Notice FSN 104121

2024-11-18

1040266224	1040271132	1040274187	1040278717
1040266428	1040271134	1040274257	1040278923
1040266437	1040271135	1040274260	1040278972
1040266595	1040271136	1040274299	1040278979
1040266596	1040271140	1040274301	1040279258
1040266597	1040271141	1040274302	1040279361
1040266598	1040271144	1040274303	1040279381
1040266600	1040271145	1040274305	1040279580
1040266760	1040271637	1040274306	1040279938
1040266886	1040271650	1040274309	1040279950
1040267279	1040271651	1040274416	1040279962
1040267280	1040271695	1040274417	1040280325
1040267794	1040271698	1040274799	1040280873
1040268108	1040271699	1040274867	1040280905
1040268342	1040271892	1040274969	1040280906
1040268909	1040271976	1040274970	1040280908
1040268947	1040271977	1040275350	1040280909
1040269173	1040271993	1040275355	1040280925
1040269174	1040272006	1040275369	1040280966
1040269294	1040272009	1040275564	1040280981
1040269404	1040272015	1040275568	1040281021
1040269459	1040272031	1040275687	1040281028
1040269775	1040272034	1040275856	1040281034
1040270073	1040272215	1040275967	1040281037
1040270078	1040272220	1040276020	1040281040
1040270242	1040272236	1040276059	1040281042
1040270276	1040272273	1040276904	1040281125
1040270277	1040272313	1040276905	1040281581
1040270437	1040272723	1040276906	1040282003
1040270442	1040272744	1040277321	1040282054
1040270604	1040273099	1040277323	1040282084
1040270605	1040273100	1040277378	1040282085
1040270668	1040273270	1040277761	1040282086
1040270806	1040273271	1040277762	1040282234
1040270835	1040273615	1040277825	1040282280
1040270879	1040273635	1040278243	1040282296
1040270902	1040273637	1040278258	1040282297
1040270903	1040273639	1040278342	1040282305
1040270966	1040273640	1040278347	1040282366
1040270982	1040273827	1040278523	1040282405
1040270983	1040273829	1040278524	1040282415
1040271019	1040273830	1040278543	1040282437
1040271022	1040273831	1040278544	1040282637
1040271050	1040274125	1040278564	1040282638
1040271057	1040274181	1040278711	1040282699
1040271123	1040274182	1040278713	1040282704

Important Safety Notice to the Market

Safety Notice FSN 104121

2024-11-18

1040282709	1040289265	1040297066	1040302903
1040282710	1040289266	1040297069	1040302904
1040282711	1040289465	1040297071	1040302905
1040282712	1040290023	1040297073	1040302906
1040283004	1040290027	1040297449	1040302932
1040283005	1040290030	1040297672	1040302933
1040283006	1040290177	1040297914	1040302964
1040283443	1040291174	1040298481	1040303179
1040283550	1040291237	1040298482	1040303182
1040283732	1040291238	1040298489	1040303419
1040283740	1040291239	1040298801	1040303420
1040283741	1040291240	1040298812	1040303421
1040284397	1040291241	1040298845	1040303422
1040284463	1040291242	1040298846	1040303423
1040284804	1040291243	1040298847	1040303424
1040284805	1040291269	1040298848	1040303428
1040284849	1040291383	1040298849	1040303932
1040284934	1040291573	1040298850	1040304034
1040284936	1040291591	1040298851	1040304036
1040284937	1040291592	1040298852	1040304038
1040285027	1040291593	1040298853	1040305635
1040285173	1040292115	1040298854	1040305636
1040285266	1040292116	1040299036	1040305657
1040285274	1040292117	1040299074	1040305845
1040285279	1040292118	1040299075	1040305860
1040285281	1040292119	1040299091	1040305862
1040285331	1040292142	1040299123	1040305917
1040285333	1040292412	1040299125	1040307280
1040285453	1040292454	1040299126	1040307281
1040286083	1040292621	1040299127	1040307282
1040286522	1040295975	1040299128	1040307283
1040286567	1040295977	1040299129	1040307300
1040286868	1040295988	1040299162	1040307399
1040286906	1040296010	1040299192	1040308068
1040286907	1040296013	1040299198	1040308072
1040287396	1040296108	1040299787	1040308081
1040288401	1040296383	1040300320	1040308082
1040288597	1040296594	1040300325	1040308083
1040288598	1040296599	1040302502	1040308084
1040288721	1040296604	1040302894	1040308085
1040289106	1040296893	1040302897	1040308086
1040289107	1040297026	1040302898	1040308140
1040289186	1040297031	1040302899	1040308471
1040289188	1040297039	1040302900	1040308996
1040289262	1040297044	1040302901	1040309031
1040289264	1040297060	1040302902	1040309032

Important Safety Notice to the Market

Safety Notice FSN 104121

2024-11-18

1040309044	1040318437	1040324108	1040330019
1040309045	1040318438	1040324109	1040330021
1040309046	1040318439	1040324110	1040330022
1040309047	1040318440	1040324111	1040330024
1040309048	1040318441	1040324117	1040330946
1040309049	1040318454	1040324951	1040331208
1040309050	1040318456	1040324955	1040331210
1040309051	1040318874	1040324956	1040332516
1040309052	1040319381	1040324957	1040332518
1040309053	1040319385	1040324958	1040332681
1040309055	1040319635	1040324959	1040332684
1040309056	1040319659	1040324960	1040332749
1040309061	1040320130	1040324963	1040332783
1040309063	1040320132	1040324966	1040332969
1040309080	1040320155	1040324969	1040332970
1040309244	1040320576	1040327008	1040334847
1040309245	1040320681	1040327009	1040334933
1040309246	1040320693	1040327233	1040334955
1040309253	1040320694	1040327301	1040335039
1040309414	1040321593	1040328162	1040335045
1040309415	1040321863	1040328261	1040335047
1040309457	1040321864	1040328262	1040336725
1040309458	1040321865	1040328303	1040336831
1040309459	1040321866	1040328315	1040336832
1040309530	1040321868	1040328419	1040336833
1040310992	1040321869	1040328598	1040336834
1040311070	1040321879	1040328599	1040336835
1040313699	1040322109	1040329334	1040336836
1040313901	1040322110	1040329432	1040340271
1040314171	1040322143	1040329433	1040343166
1040314204	1040322144	1040329434	1040343723
1040314207	1040322416	1040329493	1040343794
1040314208	1040322565	1040329518	1040345233
1040314772	1040322988	1040329539	1040345234
1040314773	1040323017	1040329543	1040345235
1040314776	1040323030	1040329557	1040346433
1040314792	1040323690	1040329558	1040348849
1040314793	1040323691	1040329672	1040352309
1040314869	1040323692	1040329676	1040352741
1040315233	1040323693	1040329677	1040352753
1040315665	1040323729	1040329906	1040354501
1040316672	1040323732	1040330011	1040355871
1040316674	1040323734	1040330012	1040355872
1040316676	1040324096	1040330013	1040356817
1040316732	1040324105	1040330014	1040362721
1040318432	1040324107	1040330015	1040363481

Important Safety Notice to the Market

Safety Notice FSN 104121
2024-11-18

1040363602

1040380545